


CATALOG NUMBER	MATERIALS	
	BOTTLE	CLOSURE
2241-0050	HDPE	PP

NOMINAL CAPACITY (L/GAL)	25/6.5	
NECK FINISH (mm)	120-415	
THREAD TPI	5	
NOMINAL BOTTLE WEIGHT (oz/g)	68/1928	
NOMINAL CLOSURE WEIGHT (oz/g)	1.4/40	
	INCH	MM (REF)
A (NECK ID)	4.25	108
B (HEIGHT W/CLOSURE)	17.76	451
C (HEIGHT W/O CLOSURE)	17.32	440
D (SHOULDER HEIGHT)	14.25	362
E (LENGTH)	12.25	311
F (WIDTH)	8.75	222
W (MIN WALL THICKNESS)	.04	1.0



SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

### CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE P. PELYCHATY	APP'D PWC	PSS 2241-0050	REV. 0816
CK'D JC	8/29/16		
DESCRIPTION: HEAVY DUTY WIDE MOUTH JUG, 25L/6.5GAL			