


CATALOG NUMBER	NOMINAL CAPACITY (mL)	MATERIALS	
		BOTTLE	CLOSURE
2423-0500	500	LDPE	PP WHITE, YELLOW, BLUE, RED (4 EACH)

SIZE NOM. CAP (mL/oz):	500/16	
NOM. BRIM CAP (mL)...	580	
NECK FINISH (mm).....:	38-415	
THREAD (T.P.I.).....:	6	
NOM. WEIGHT WITH CLOSURE (grams).....:	58	
	INCH	MM
A (NECK I.D.).....:	1.16 ±0.03	29.5 ±0.8
B (HT. W/ CLOSURE).....:	6.81 ±0.09	173.0 ±2.3
C (HT. W/O CLOSURE)....:	6.69 ±0.09	170.0 ±2.3
D (DIA.).....:	3.02 ±0.08	76.7 ±2.0
W.T. (NOM. WALL THK.):	0.040	1.0



SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	TL	APP'D	GRP	PSS 2423-0500	REV. 0310
CK'D	PB	3/11/10			
DESCRIPTION: Color-Coded Unitary Wash Bottles					