



CATALOG NUMBER.....	SEE BELOW	
MATERIAL.....	PS	
	INCH	MM
A (BASE LENGTH).....	5.04	128.0
B (BASE WIDTH).....	3.38	85.7
C (OVERALL HEIGHT).....	.22	5.6
D (SIDEWALL HEIGHT).....	.18	4.6
E (CUTOUT, TOP; LENGTH)....	3.30	83.7
F (CUTOUT, BOTTOM; LENGTH)..	3.46	87.8
G (CUTOUT, TOP; WIDTH).....	2.17	55.2
H (CUTOUT, BOTTOM; WIDTH)..	2.33	59.2
I (CUTOUT, HEIGHT).....	0.08	2.0

CATALOG NUMBER	DESCRIPTION
253623	NON-STERILE
253624	STERILE

FIXED - 253624 WAS 263624  
J00 - 12/29/04

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE MC <sup>2</sup>	APP'D PVC	PSS 25362X-2	REV. 1002
CHK'D RDR 11/1/02	10/31/02		

DESCRIPTION:

LOW BOT UNIVERSAL LID - NUNC