



CATALOG NUMBER.....	264623	
MATERIAL.....	PS - CLEAR	
	INCH	MM
A (LENGTH).....	5.02	127.5
B (WIDTH).....	3.36	85.3
C (HEIGHT W/O STACKING LEGS):	0.31	7.9
D (OVERALL HEIGHT).....	0.35	9.0



TOP VIEW



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	JO	APP'D	BR	PSS 264623	REV.
CK'D	AEG		6/14/04		0604
DESCRIPTION:					
96 WELL LID - STERILE					