


| | | |
|----------------------------|------------|------|
| CATALOG NUMBER | SEE BELOW | |
| WORKING RANGE (µL) | 5 – 240 µL | |
| MATERIAL | PP | |
| | MM | INCH |
| A, BASE LENGTH | 127,8 | 5,03 |
| B, BASE WIDTH | 85,5 | 3,37 |
| C, TOP LENGTH | 123,5 | 4,87 |
| D, TOP WIDTH | 81,2 | 3,20 |
| E, FLANGE LONG SIDE | 2,5 | 0,10 |
| F, OVERALL HEIGHT | 22 | 0,87 |
| G1, TOP | 3,9 | 0,15 |
| G2, BOTTOM | 3,4 | 0,13 |
| H, A - I LOCATION (X) | 12,1 | 0,48 |
| J, A - I LOCATION (Y) | 9 | 0,35 |
| K, WELL TO WELL SPACING | 4,5 | 0,18 |
| L, WELL DEPTH | 19,5 | 0,77 |
| M, BOTTOM OF WELL DISTANCE | 2,5 | 0,10 |
| N, FLANGE SHORT SIDE | 2,5 | 0,10 |

| CATALOG NO | SURFACE | COLOR | STERILE | WITH LID* | UNITS/ PK/CASE |
|------------|-------------|---------|---------|-----------|----------------|
| 269390 | NON-TREATED | NATURAL | - | - | 5/60 |
| | | | | | |

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATION PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



ThermoFisher
SCIENTIFIC

Roskilde Site · Kamstrupvej 90 · P.O.Box 280
DK-4000 Roskilde · Denmark
Tel +45 4631 2000 · Fax +45 4631 2042
E-mail: info.nunc@thermofisher.com

www.nuncbrand.com

SUBJECT:

384 DEEP WELL PLATE, PP

DRAWN BY/DATE:

CHANGE ORDER NO.:

VALID FROM:

APPROVED BY/DATE:

FORM

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC