


MATERIAL	PC	
SIZE (NOMINAL VOLUME mL)	15	
NOMINAL BRIM CAPACITY (mL)	16	
NOMINAL WEIGHT (g)	7	
	IN	MM (REF)
A (TOP INSIDE DIAMETER)	.59±.01	15.0
B (BODY DIAMTER)	.66±.01	16.8
C (HEIGHT)	4.72±.02	119.9

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



GRADUATIONS ARE IN
0.1mL INCREMENTS

SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE P. PELYCHATY	APP'D 4/28/16
CK'D JC	PWC

PSS 3105-0015

REV.
0416

DESCRIPTION:

CONICAL CENTRIFUGE TUBE, 15mL