


MATERIAL	PC	
SIZE (NOMINAL VOLUME mL)	50	
NOMINAL BRIM CAPACITY (mL)	55	
NOMINAL WEIGHT (g)	15	
	IN	MM (REF)
A (TOP INSIDE DIAMETER)	1.05±.01	26.7
B (BODY DIAMTER)	1.13±.01	28.7
C (HEIGHT)	5.29±.02	134.4



SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

**CONFIDENTIAL**

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE P. PELYCHATY	APP'D 4/28/16	PSS 3105-0050	REV. 0416
CK'D JC	PWC		
DESCRIPTION: CONICAL CENTRIFUGE TUBE, 50mL			