


BOTTLE MATERIAL	AMBER HDPE	
CLOSURE MATERIAL	AMBER PP	
SIZE (NOMINAL VOLUME)	60 ML	
NOMINAL BRIM CAPACITY	69 ML	
NECK FINISH	20MM	
THREAD TPI	8	
NOMINAL WEIGHT BOTTLE W/CLOSURE	12 GRAMS	
	INCH	MM
A (CLOSURE OD)	.96	24.4
B (NECK OD)	.70	17.6
C (NECK ID)	.54	13.7
D (DIMETER)	1.54	39.1
E (HEIGHT WITHOUT CLOSURE)	3.24	82.3
F (HEIGHT WITH CLOSURE)	3.62	91.9
G (NOMINAL WALL THICKNESS)	.025	.64



SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE ED LEE	APP'D PWC	PSS 312004-0002	REV. 0317
CK'D JC			
DESCRIPTION: 60ml, NARROW MOUTH BOTTLE			