

CATALOG NO.	NOMINAL CAPACITY, (oz)	MATERIALS	
		BOTTLE	CLOSURE
312088-0016	16	LDPE	PP


SIZE CODE.....	-0016	
SIZE (NOM VOL mL/oz).....	500 / 16	
NOM. BRIM CAP. (mL).....	555	
NECK FINISH (mm).....	28-415	
THREAD (T.P.I.).....	6	
NOM. WEIGHT W/O CLOSURE (gms).....	40	
NOM. WEIGHT W/ CLOSURE (gms).....	44	
	INCH	MM
T (MAJOR THREAD DIA).....	1.09 +.00/-.03	27.7 +0.0/-0.76
E (MINOR THREAD DIA).....	0.99 +.00/-.03	25.1 +0.0/-0.76
I (NECK I.D.).....	0.84 ±.02	21.3 ±0.5
S (THREAD START).....	0.09 ±.02	2.3 ±0.5
H (FINISH HEIGHT).....	0.87	22.1
L (BEAD LOCATING DIMENSION).....	0.56	14.2
W (BEAD WIDTH).....	0.09	2.3
B (BEAD DIAMETER).....	1.28	32.5
C (HT. W/O CLOSURE).....	6.62 ±.09	168.2 ±2.3
Z (HT. W/ CLOSURE).....	6.78 ±.09	172.2 ±2.3
D (DIAMETER OF BODY).....	2.88 ±.08	73.2 ±2.0
W.T. (NOMINAL WALL THK).....	0.04	1.02
MINIMUM WALL THK.....	0.015	0.38

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

SPECS ON BACK

NOT TO SCALE



CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	KOW	APP'D	GRP	PSS 31-2088-0016	REV.
CK'D	TL		4/10/12		0412
DESCRIPTION:				500ml NARROW MOUTH PACKAGING BOTTLE	