


CATALOG NO.	NOMINAL CAPACITY, (oz)	MATERIALS	
		BOTTLE	CLOSURE
312084-0004	4	TRANSLUCENT AMBER HOPE	OPAQUE AMBER PP
312087-0004	4	PPCO	PP

SIZE CODE	-0004	
SIZE (NOM VOL mL/oz)	125/4	
NOM. BRIM CAP. (mL):	143	
NECK FINISH (mm)....	24-415	
THREAD (T.P.I.).....	8	
NOM. WEIGHT W/CLOSURE (gms).....	16	
	INCH	MM
A (NECK I.D.).....	0.69 ± 0.02	17.5 ± 0.5
B (HT. W/CLOSURE)...	3.97 ± 0.08	100.8 ± 2.0
C (HT. W/O CLOSURE):	3.87 ± 0.08	98.3 ± 2.0
D (DIA).....	1.99 ± 0.06	50.5 ± 1.5
W.T.(NOM. WALL THK):	0.030	0.76
T (DIA).....	.930 ± 0.020	23.62 ± 0.51
E (DIA).....	.848 ± 0.015	21.50 ± 0.40
H.....	.790 ± 0.020	20.00 ± 0.50
I (DIA).....	.735 ± 0.020	18.70 ± 0.50



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	TL	APP'D	PSS 31-208X-0004-S	REV.
CK'D	DM	GRP		1209

DESCRIPTION:

NARROW-MOUTH SAMPLE BOTTLES