


CATALOG NO.....:	312185-0008
NOM. CAP. (oz).....:	8
MATERIALS	BOTTLE...: OPAQUE AMBER HDPE
	CLOSURE...: OPAQUE AMBER PP
SIZE (NOM VOL mL/oz)	250/8
NOM. BRIM CAP. (mL):	300
NECK FINISH (mm)....:	43-415
THREAD (T.P.I.).....:	6
NOM. WEIGHT W/CLOSURE (gms).....:	35
A (NECK I.D.).....:	1.300 ±.020
B (HT. W/CLOSURE)...:	5.200 ±.045
C (HT. W/O CLOSURE):	5.025 ±.045
D (DIA).....:	2.400 ±.045
W.T.(NOM WALL THK).:	.025

THIS DRAWING IS NOT TO BE REPRODUCED  
-CONFIDENTIAL-

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MCC	APP'D	PVC	PSS 31-2185-0008-S	REV.	-
CK'D	PVC		3/6/92			
DESCRIPTION:				WIDE-MOUTH, ROUND PACKAGING BOTTLES		