


Technical Data Sheet

Connect Edge: 3.7 Meter Door Switch Sensor

Revision A

Specifications	Catalog Number and Production Model Number	EDGEPTSENDU004	
Application, Rating and Electrical Data			
Application	Sensor used for monitoring door open/close		
Component Type	Cable: PVC; PVC shrink tubing <ul style="list-style-type: none"> Switch Dimensions: 13 mm x 14mm x 100 mm Magnet Dimensions: 13 mm x 14mm x 100 mm ABS Plastic 		
Module and Sensor Information			
Operation Gap	70mm		
Cable Length with connector	3.7 meters		
Switch cycles	50 million		
Weight	160g		
General			
Mounting Options	Connect Edge Adapter Kit, 6-Sensor Configuration, Analog Ports		
Disinfection/Cleaning	<ul style="list-style-type: none"> Clean the Sensor with a soft cloth lightly moistened with water, a detergent or isopropanol alcohol. Do not use any aggressive cleaning agents or scratching cleansers that might damage your sensor 		
Memory			
Sampling Rates	Once per minute		
Standards and Regulations	ROHS, FCC, CE, TUV-SUD		
Environmental Conditions			
Sensor Operating Conditions	-25°C to 70°C (-15°F to 160°F)		
Country of Origin	China		
Warranty			
Sensor Warranty	1 year		