

The Thermo Scientific Barnstead Accudispense provides the ultimate method for dispensing high purity water from your Thermo Scientific Barnstead Nanopure® or TII water systems.

Thermo Scientific Barnstead Accudispense™



Easy dispensing up to eight feet (2.5 meters) away



Includes 0.2 micron absolute final filter



Product Description

Dispense water remotely with our patented volumetric dispense feature. Allows you to mount the Nanopure or TII virtually anywhere in your lab while maintaining a source of high purity water and convenient dispensing at a location you choose.

Product Features

- Accurately dispense between 250 ml and 60 L into a flask, cylinder or carboy.
- Effortlessly dispense water up to eight feet (2.5 meters) away from the Nanopure or TII.

- LED on the base indicates that the water is above or below your programmed purity level.
- Communicates directly with the Nanopure or TII allowing you to initiate the dispense from the bench.
- Includes a 0.2 micron absolute final filter to remove bacteria and particles while maintaining flow.
- Recirculating up to the tip of the dispenser eliminates contamination of your water.



LED display for programmed water purity level



Connects directly to your Nanopure

Thermo Scientific Barnstead Accudispense

Technical Specifications		
Model Number	Description	Shipping Weight
D13661	Accudispense Volumetric Remote Dispenser	37 lbs. (16.78 kg)
D3750	0.2 Micron Final Filter	1 lb. (0.45 kg)

North America: USA/Canada +1 800 553 0039

Europe: Austria +43 1 801 40 0, Belgium +32 2 482 30 30, France +33 2 2803 2000, Germany national toll free 08001-536 376, Germany international +49 6184 90 6940, Italy +39 02 02 95059 341, Netherlands +31 76 571 4440, Nordic countries +358 9 329 100, Russia/CIS +7 (812) 703 42 15, Spain/Portugal +34 93 223 09 18, Switzerland +41 44 454 12 12, UK/Ireland +44 870 609 9203

Asia: China +86 21 6865 4588 or +86 10 8419 3588, India +91 22 6716 2200, Japan +81 45 453 9220, Other Asian countries +852 2885 4613

Countries not listed: +49 6184 90 6940 or +33 2 2803 2000

www.thermo.com/purewater

© 2008 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

DS-LEWP-ACUDIPSE-0608

Thermo
SCIENTIFIC