



# Improve supply chain availability of sterile products faster

## Switch from gamma to e-beam sterilization of your Thermo Scientific™ Nalgene™ products and realize the benefits:

- **Faster turnaround**—less time required to process your sterile products
- Identical sterility assurance levels (SAL)
- Identical validation standard (ISO 11137)
- RNase- and DNase-free and nonpyrogenic, as always; certified DNA-free where applicable
- Does not require a dangerous radioactive isotope like gamma; e-beam is an on/off electricity-based technology with a faster sterilization process than gamma



## Ordering information

Product	Volume	Unit size	Cat. No.
Nalgene PPCO Micro Packaging Vials with E-beam Irradiation	0.5 mL	Case of 1,000	342810-0005
	2 mL	Case of 1,000	342810-0020

Find out more at [thermofisher.com/ebeam-MPV](https://thermofisher.com/ebeam-MPV)