




CATALOG NO.	NOMINAL TUBE SIZE, (ml)		MATERIAL
005973	0.5	1.5	RESMER

TUBE SIZE.....:	0.5		1.5	
ARRAY.....:	4 x 6		4 x 6	
	INCHES	mm	INCHES	mm
A (LENGTH).....:	5.73±.16	145.5±4.0	5.73±.16	145.5±4.0
B (WIDTH).....:	3.62±.14	91.9±3.6	3.64±.14	92.4±3.6
C (HEIGHT).....:	2.15±.10	54.6±2.5	2.15±.10	54.6±2.5
D (DIA.).....:	0.33±.01	8.4±0.3	0.44±.01	11.2±0.3
NOM.WT.(gm).....:	85		80	
COLOR	SIZE CODE			
WHITE	-0005		-0015	

TUBE SIZE.....:	0.5		1.5	
ARRAY.....:	8 × 12		8 × 12	
	INCHES	mm	INCHES	mm
A (LENGTH).....:	6.75±.16	171.5±4.0	7.90±.18	200.7±4.6
B (WIDTH).....:	4.42±.14	112.3±3.6	5.18±.15	131.6±3.8
C (HEIGHT).....:	2.15±.10	54.6±2.5	2.15±.10	54.6±2.5
D (DIA.).....:	0.33±.01	8.4±0.3	0.44±.01	11.2±0.3
NOM.WT.(gn).....:	105		140	
COLOR	SIZE CODE			
WHITE	-9005		-9015	
BLUE	-9305		-9315	
RED	—		-9515	



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE IN CHINA	MC	APP'D PVC 1/31/03	PSS 00-5973-XXXX-1	REV. 0103
DESCRIPTION: MICROCENTRIFUGE TUBE RACKS				