


NOMINAL CAPACITY	500mL	
BRIM CAPACITY	631mL	
NOMINAL WEIGHT	108g	
	INCH	MM (REF)
A (NOMINAL DIAMETER)	2.05	52.1
B (HEIGHT)	14.25	362.0
C (BASE DIAMETER)	4.63	117.6
D (NOMINAL WALL THICKNESS TOP)	.065	1.7
E (NOMINAL WALL THICKNESS BOTTOM)	.085	2.2

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	PMP	APP'D	9/26/16	PSS 6230-0500	REV. 0916
CK'D	JC	PWC			
DESCRIPTION: 500mL HYDROMETER JAR					