


CATALOG NUMBER.....:	4000-0050	4000-0100	4000-0200			
NOM. CAP. (mL).....:	50	100	200			
CLOSURE SIZE (mm)....:	13	20	20			
NOMINAL WT WITH CLOSURE (gms).....:	14	22	29			
MATERIAL	FLASK.....:	PP	PP			
	CLOSURE...:	PP	PP			
	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D. AT TOP.):	.34 ±.04	8.6 ±1.0	0.52 ±.04	13.2 ±1.0	0.52 ±.04	13.2 ±1.0
B (NECK O.D.).....:	.47 ±.04	11.9 ±1.0	0.59 ±.04	15.0 ±1.0	0.66 ±.04	16.8 ±1.0
D (BASE DIA).....:	1.38 ±.04	35.1 ±1.0	1.21 ±.04	30.7 ±1.0	2.10 ±.04	53.3 ±1.0
C (BULB DIA).....:	1.82 ±.06	46.2 ±1.5	2.18 ±.06	55.4 ±1.5	2.75 ±.06	69.9 ±1.5
E (HT W/O CLOSURE)...:	5.69 ±.12	144.5 ±3.0	7.34 ±.12	186.4 ±3.0	8.25 ±.12	209.5 ±3.0
F (HT WITH CLOSURE)..:	5.8 ±.2	147 ±5	7.5 ±.2	191 ±5	8.4 ±.2	213 ±5

CATALOG NUMBER.....:	4000-0250	4000-0500	4000-1000			
NOM. CAP. (mL).....:	250	500	1000			
CLOSURE SIZE (mm)....:	20	24	28			
NOMINAL WT WITH CLOSURE (gms).....:	41	84	123			
MATERIAL	FLASK.....:	PP	PP			
	CLOSURE...:	PP	PP			
	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D. AT TOP.):	0.52 ±.04	13.2 ±1.0	0.70 ±.04	17.8 ±1.0	0.86 ±.04	21.8 ±1.0
B (NECK O.D.).....:	0.72 ±.04	18.3 ±1.0	0.81 ±.04	20.6 ±1.0	0.94 ±.04	23.9 ±1.0
C (BASE DIA).....:	2.25 ±.04	57.2 ±1.0	2.75 ±.04	69.9 ±1.0	3.62 ±.04	91.9 ±1.0
D (BULB DIA).....:	3.00 ±.06	76.2 ±1.5	3.75 ±.06	95.3 ±1.5	4.75 ±.06	120.7 ±1.5
E (HT W/O CLOSURE)...:	8.56 ±.12	217.4 ±3.0	10.44 ±.12	265.2 ±3.0	12.19 ±.12	309.6 ±3.0
F (HT WITH CLOSURE)..:	8.9 ±.2	226 ±5	10.8 ±.3	274 ±8	12.4 ±.3	315 ±8

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	PMP	APP'D	PWC	PSS 4000-XXXX	REV.
CK'D	JC		9/7/16		0916

DESCRIPTION:

VOLUMETRIC FLASKS