



Loading and compatibility guide for Banded Analytical Membranes

Thermo Scientific™ Banded Analytical Membranes for Dispensers are compatible with multiple hardware platforms for loading in one of two possible orientations. Table 1 references common dispenser hardware and the direction the box should be facing to properly load the membranes.* Note the orientation arrows on each membrane box.

Table 1. Compatible hardware platforms according to box orientation.

Compatible hardware for dispensing from the side	Compatible hardware for dispensing from the top
	
<p>MilliporeSigma™ EZ Pak™ and EZ Pak™ Curve membrane dispensers Note: Other hardware platforms may be compatible.</p>	<p>Sartorius™ Microsart™ and Microsart™ e.motion dispensers Note: Other hardware platforms may be compatible.</p>

* Follow the standard installation instructions for your hardware platform.

Learn more at thermofisher.com/filtration