



Specifications	Model Number				Representative Model (tested)
	TSX1205PA, TSX1205PD, TSX1205PY, TSX1205PZ, CH1205PA, CH1205PD				
	Recommended cabinet temperature setpoint		Capacity m ³		
	5°C		0.19		
	Application				
Vaccine storage					
Program	Temperature requirements				Representative Model (tested)
	Steady state	Short door openings	Long (3 min.) door opening recovery		
NSF 456	2°C < T < 8°C	2°C < T < 8°C	in < 15min, < 8°C		TSX1205PA

