

Specifications	Model Number				Representative Model (tested)
	TSX2305PA, TSX2305PD, TSX2305PY, TSX2305PZ, CH2305PA, CH2305PD				
	Recommended cabinet temperature setpoint		Capacity m ³		
	5°C		0.36		
	Application				
Vaccine storage					
Program	Temperature requirements				Representative Model (tested)
	Steady state	Short door openings	Long (3 min.) door opening recovery		
NSF 456	2°C < T < 8°C	2°C < T < 8°C	in < 15min, < 8°C		TSX2305PA

