



Specifications	Model Number				Representative Model (tested)
	TSX3005PA, TSX3005PD, TSX3005PY, TSX3005PZ, CH3005PA, CH3005PD				
	Recommended cabinet temperature setpoint		Capacity m <sup>3</sup>		
	5°C		0.46		
	Application				
Vaccine storage					
Program	Temperature requirements				Representative Model (tested)
	Steady state	Short door openings	Long (3 min.) door opening recovery		
NSF 456	2°C < T < 8°C	2°C < T < 8°C	in < 15min, < 8°C		TSX3005PA

