

Specifications	Model Number				
	TSX3030FA, TSX3030FD, TSX3030FY, TSX3030FZ, CH3030FA, CH3030FD				
	Recommended cabinet temperature setpoint			Capacity m <sup>3</sup>	
	-30°C			0.59	
	Application Vaccine storage				
Program	Temperature requirements				Representative Model (tested)
	Steady state	Short door openings	Long (1 min.) door opening recovery	Defrost	
NSF 456	T < -15°C	T < -15°C	in < 15min, < -15°C	in < 15min, < -15°C	TSX3030FA

