

Specifications	Model Number				Representative Model (tested)
	TSX5005PA, TSX5005PD, TSX5005PY, TSX5005PZ, CH5005PA, CH5005PD				
	Recommended cabinet temperature setpoint		Capacity m ³		
	5°C		0.71		
	Application				
Vaccine storage					
Program	Temperature requirements				
	Steady state	Short door openings	Long (3 min.) door opening recovery		
NSF 456	2°C < T < 8°C	2°C < T < 8°C	in < 15min, < 8°C	TSX5005PA	

