

Specifications	Model Number				
	TSX5030FA, TSX5030FD, TSX5030FY, TSX5030FZ, CH5030FA, CH5030FD				
	Recommended cabinet temperature setpoint			Capacity m ³	
	-30°C			0.92	
	Application Vaccine storage				
Program	Temperature requirements				Representative Model (tested)
	Steady state	Short door openings	Long (1 min.) door opening recovery	Defrost	
NSF 456	T < -15°C	T < -15°C	in < 15min, < -15°C	in < 15min, < -15°C	TSX5030FA

