




|                               |           |        |
|-------------------------------|-----------|--------|
| CATALOG NUMBER.....:          | SEE BELOW |        |
| NOM. CAP. (mL).....:          | 0.35      |        |
| MATERIAL.....:                | PP        |        |
|                               | INCH      | MM     |
| A (BASE LENGTH).....:         | 5.030     | 127.76 |
| B (BASE WIDTH).....:          | 3.365     | 85.47  |
| C (TOP LENGTH).....:          | 4.85      | 123.2  |
| D (TOP WIDTH).....:           | 3.18      | 80.8   |
| F (OVERALL HEIGHT).....:      | 0.57      | 14.5   |
| G DIA.....:                   | 0.14      | 3.6    |
| H (A-1 LOCATION) (X).....:    | 0.47      | 11.9   |
| J (A-1 LOCATION) (Y).....:    | 0.35      | 8.9    |
| K (WELL-TO-WELL SPACING)....: | 0.18      | 4.6    |
| L (WELL DEPTH).....:          | 0.26      | 6.6    |
| N (FLANGE SHORT SIDE).....:   | 0.10      | 2.5    |

|                |                   |
|----------------|-------------------|
| CATALOG NUMBER | COLOR             |
| 267459         | NATURAL           |
| 267460         | NATURAL (STERILE) |
| 267461         | BLACK             |
| 267462         | WHITE             |



TOP VIEW



SECTION X-X

SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

|                     |    |              |         |            |              |
|---------------------|----|--------------|---------|------------|--------------|
| MADE<br>IN<br>CHINA | PB | APP'D<br>PVC | 1/18/05 | PSS 2674XX | REV.<br>0105 |
|---------------------|----|--------------|---------|------------|--------------|

DESCRIPTION:  
384, 35µl, PP, WELL PLATE

- NOTES:
1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
  2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.