


CATALOG NUMBER	SEE BELOW	
WORKING RANGE. (µL)	50 – 250 µL/WELL	
MATERIAL	PS	
	MM	INCH
A , BASE LENGTH	127,8	5,03
B , BASE WIDTH	85,6	3,37
C , TOP LENGTH	123,5	4,86
D , TOP WIDTH	82,3	3,2
E , FLANGE LONG SIDE	2,7	0,11
F , OVERALL HEIGHT	14,4	0,60
G1 , DIA TOP	Ø7	Ø0,28
G2 , DIA BOTTOM	Ø6,3	Ø0,25
H , A - I LOCATION (X)	14,5	0,57
J , A - I LOCATION (Y)	11,5	0,45
K , WELL TO WELL SPACING	9	0,35
L , WELL DEPTH	10,3	0,41
M , BOTTOM OF WELL DISTANCE	4,1	0,16
N , FLANGE SHORT SIDE	2,7	0,11

CATALOG NO	SURFACE	COLOR	STERILE	WITH LID*	UNITS/PK/CASE
280210	NON-TREATED	CLEAR	+	-	10/180

NOTES:

- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
- THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATION PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



ThermoFisher
SCIENTIFIC

Roskilde Site • Kamstrupvej 90 • P.O.Box 280
DK-4000 Roskilde • Denmark
Tel +45 4631 2000 • Fax +45 4631 2042
E-mail: info.nunc@thermofisher.com

www.nuncbrand.com

SUBJECT:

U96 MICROWELL PLATE, PS

DRAWN BY/DATE:

CHANGE ORDER NO.:

VALID FROM:

APPROVED BY/DATE:

FORMA

ITEM NUMBER:

MATERIAL:

DRAWING NO.:

VERSION:

PAGE TYPE:

SCALE:

METRI

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC