

# Thermo Scientific Tube Revolver

Increase your lab's efficiency with our variable speed tube revolver. Its small footprint, flexible design, and simple display help increase flexibility and effectiveness.

## FEATURES

- Interchangeable rotisseries accommodate different tube sizes for maximum flexibility. Each paddle can be adjusted from 0-90 degrees.
- Achieve more effective mixing with "oscillating" mode. Oscillation will occur at every 40 degrees in the rotation cycle.
- Small footprint
- Plug-and-play simplicity
- Maintenance-free motor, quiet operation



## APPLICATIONS

Molecular Biology | Clinical  
Biochemistry | Others

4

### Ordering Information and Specifications

Cat. No.	88881001	88881002
Plug Configurations		
Rotation Speed	10-40 rpm	
Capacity	<b>Combination of 4 paddles of:</b> 7 x 10/15ml tubes 7 x 5/7ml tubes 15 x 1.5ml/2ml tubes 21 x 0.5/0.8ml tubes One paddle of 6 x 50ml tubes	
Standard Rotisserie Paddles included	88881005 – 10/15ml paddle, quantity 2 88881006 – 5ml/7ml paddle, quantity 2 88881007 – 1.5ml/2ml paddle, quantity 2 88881009 – 50ml paddle, quantity 1  Paddles can also be purchased as accessories in quantities listed above	
Operating Temperature	4°C – 60°C	
Dimensions W x L x H	225 x 154 x 195mm / 8.9 x 6 x 7.7in	
Max Power	6W	
Weight	1.2kg / 2.65lbs	
Power	AC100V-240V / 50Hz – 60Hz	
Certification	cCSAus, CE	
Warranty	2 years	

### Ordering Information and Specifications

	Cat. No.
Paddles 10ml/15ml, 2 pack	88881005
Paddles 5ml/7ml, 2 pack	88881006
Paddles 1.5ml/2ml, 2 pack	88881007
Paddles 0.5ml/0.8ml, 2 pack	88881008
Paddle 50ml	88881009



Paddles